# Sil-Select Plus

# Sil-Select Plus

## Sil-Select Plus™

Gradient system for semen preparation

STERILE A

Sterilized by sterile filtration. Document reference: FP09 I13 R01 C.1 Update: 05.02.2018

### USED ABBREVIATIONS

ICSI Intracytoplasmatic Sperm Injection  
IVF In Vitro Fertilization  
IUI Intra Uterine Insemination

### INTENDED USE AND COMPOSITION

Sil-Select Plus is a ready-to-use gradient system for semen preparation. The following gradients are available: – Sil-Select Plus Upper layer (45%) – Sil-Select Plus 80% – Sil-Select Plus Lower layer (90%) Sil-Select Plus gradients consist of silane-coated colloidal silica particles suspended in HEPES-buffered EBSS (Earle's balanced salt solution) supplemented with 0.4-2.2 g/l human serum albumin. The Sil-Select Plus system can be used in combination with IUI, IVF and ICSI.

For professional use only.

### MATERIAL INCLUDED WITH THE KIT

- Product code: SIP008 (8 preparation kit)  
 » 8x 2.5ml Sil-Select Plus LOWER LAYER  
 » 8x 2.5ml Sil-Select Plus UPPER LAYER  
 » 8x 5ml Sil-Select Plus Sperm Washing/Insemination medium
- Product code: SIP016 (16 preparation kit)  
 » 2x 20ml Sil-Select Plus LOWER LAYER  
 » 2x 20ml Sil-Select Plus UPPER LAYER  
 » 2x 20ml Sil-Select Plus Sperm Washing/Insemination medium
- Product code: SIP020 (2x 20ml)  
 » 1x 20ml Sil-Select Plus LOWER LAYER  
 » 1x 20ml Sil-Select Plus UPPER LAYER
- Product code: SIP050 (2x 50ml)  
 » 1x 50ml Sil-Select Plus LOWER LAYER  
 » 1x 50ml Sil-Select Plus UPPER LAYER
- Product code: SIP100 (2x 100ml)  
 » 1x 100ml Sil-Select Plus LOWER LAYER  
 » 1x 100ml Sil-Select Plus UPPER LAYER
- Product code: SIP80\_100  
 » 1x 100ml Sil-Select Plus 80%
- Product code: SIP80\_500  
 » 1x 500ml Sil-Select Plus 80%

Single bottles are available CE marked (Notified Body number 0344). Single bottles of Sil-Select Plus can be ordered separately.

### MATERIAL NOT INCLUDED WITH THE KIT

- » 3cc syringes with 1 1/2" 21 g needle
  - » Centrifuge (must be able to operate for up to 30 minutes at 400 g)
  - » Incubator or water bath at 37°C (optional)
  - » SIPO20/SIPO50/SIP100/SIP80\_100/SIP80\_500: Sil-Select Plus Sperm Washing/Insemination medium » LAF bench (ISO Class 5)
- » pH: 7.20-7.90 (Release criteria: 7.20-7.60)  
 » Osmolality: 280-310 mOsm/kg (Upper Layer) 290-330 mOsm/kg (Lower Layer)  
 » Density: 1.1050-1.1150 g/ml (Lower Layer) 290-330 mOsm/kg (80%)  
 » Viscosity: < 1.65 cP (Lower Layer, 80%)  
 » Endotoxin: < 0.5 EU/ml  
 » Sterility: SAL 10<sup>-3</sup>  
 » Sperm Survival test: ≥ 80% survival after 4 hours exposure of densely selected spermatozoa to the test medium; ≥ 75% survival after 24 hours exposure of density selected spermatozoa to the test medium  
 » Not MEA tested  
 » A certificate of analysis and MSDS are available upon request.

### CALCULATIONS OF G-FORCES

The g-force of your centrifuge can be calculated using this formula:

$$g = 1.118 \times r \times rpm^2 \div rpm = \text{Square root } (g / (1.118 \times r))$$

r = radius of centrifuge in mm  
 rpm = rotations per minute / 1000

Example 1 r = 100 mm  
 rpm = 1800 rotations per minute  
 g = 1.118 x 100 x 3.24 = 362g

Example 2 r = 100 mm  
 g = 350g  
 rpm = SQR (350 / (1.118 x 100)) = 177  
 rpm = 1770 rotations per minute

### PRE-USE CHECKS

- » Do not use the product if seal of the container is opened or defect when the product is delivered.
- » Do not use the product if it has changed colour or shows any evidence of microbial contamination.
- » Mix the density gradient bottles by 5 bottle inversions before use.

## INSTRUCTIONS FOR USE WITH FRESH SEMEN SAMPLES

- Bring all components of the system and samples to room temperature or to 37°C.
- Transfer 2.5 ml of Sil-Select Plus Upper Layer into a sterile disposable centrifuge tube.
- Using a 3cc syringe with a 1 1/2" 21 g needle, place 2.5 ml of Sil-Select Plus Lower Layer under the Upper Layer. Take care that the Upper Layer is distinctly separated. This is done by placing the tip of the needle on the bottom of the test tube and slowly dispensing the Sil-Select Plus Lower Layer. This two layer gradient is stable for up to two hours.
- Gently place up to 2.5 ml of liquefied semen onto the Upper Layer using a transfer pipette or syringe.
- Centrifuge for 15 to 18 minutes at 350 to 400 g. When this centrifugation is completed you may not be able to visibly see a pellet. If so, it is essential to continue the procedure with a second centrifugation of 3 to 5 minutes.
- Remove supernatant down to the pellet.
- Using a syringe, add 2-3 ml of sperm washing medium and re-suspend the pellet.
- Centrifuge for 8 to 10 minutes at 300 g. Higher sperm concentration will require the maximum 10 minutes centrifugation to ensure a complete and thorough sperm wash.
- Remove supernatant down to the pellet and repeat steps 7 and 8.
- Remove supernatant and replace with a suitable volume of appropriate medium.

## INSTRUCTIONS FOR USE WITH FROZEN SEMEN SAMPLES

- Bring all components of the kit and samples to room temperature or to 37°C.
- Transfer 1 ml of Sil-Select Plus Upper Layer into a sterile disposable centrifuge tube.
- Using a 3cc syringe with a 1 1/2" 21 g needle, place 1 ml of Sil-Select Plus Lower Layer under the Upper Layer. This care that the two layers are distinctly separated. This is done by placing the tip of the needle on the bottom of the test tube and slowly dispensing the Sil-Select Plus Lower Layer. This two-layer gradient is stable for up to two hours.
- Gently place a thawed semen sample onto the Upper Layer using a transfer pipette or syringe (0.5 ml maximum).
- Centrifuge for 15-20 minutes at 350 g.
- Remove supernatant down to no less than the 0.5 ml mark above the pellet.
- Using a syringe, add 2-3 ml of sperm washing medium and resuspend the pellet.
- Centrifuge for 8 to 10 minutes at 300 g.
- Remove supernatant down to the pellet and repeat steps 7 and 8.
- Remove supernatant and replace with a suitable volume of appropriate medium.

If samples do not liquefy and therefore do not pass through the layers, increasing the centrifugal force up to, but no more than, 500 g will help to separate the sperm.

## STORAGE AND CONSERVATION

- » Store products between 2-25°C before first use, once opened at 2-8°C.
- » Open and close bottles under aseptic conditions (e.g. LAF bench, ISO Class 5).
- » The products can be used safely up to 7 days after opening when sterile conditions are maintained and the products are stored at 2-8°C. Do not freeze.
- » Do not use after expiry date
- » Content cannot be re-sterilized after opening
- » The products are stable after transport (max. 5 days) at elevated temperatures (≤ 37°C)

## WARNINGS AND PRECAUTIONS

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infectious agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven viruses transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens. Always work under strict hygienic conditions (e.g. LAF-bench ISO Class 5) to avoid possible contamination.

## TECHNICAL SUPPORT

FertPro N.V.  
 Industriepark Noord 32  
 8730 Beernem / Belgium  
 Tel +32 (0)50 79 18 05  
 Fax +32 (0)50 79 17 99  
 URL: www.fertipro.com  
 E-mail: info@fertipro.com

Example 1 r = 100 mm  
 rpm = 1800 rotations per minute  
 g = 1.118 x 100 x 3.24 = 362g

Example 2 r = 100 mm  
 g = 350g  
 rpm = SQR (350 / (1.118 x 100)) = 177  
 rpm = 1770 rotations per minute

## PRE-USE CHECKS

- » Do not use the product if seal of the container is opened or defect when the product is delivered.
- » Do not use the product if it has changed colour or shows any evidence of microbial contamination.
- » Mix the density gradient bottles by 5 bottle inversions before use.

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# Sil-Select Stock

## Sil-Select Stock™

Stock solution for semen preparation

STERILE A

Sterilized by sterile filtration.  
Document reference: FP09 113 R01 C.1  
Update: 05.02.2018

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### INTENDED USE AND COMPOSITION

Sil-Select Stock is a ready-to-use isotonic gradient for semen preparation with a density of approximately 1.12 g/ml. Sil-Select Stock consists of silane-coated colloidal silica particles suspended in HEPES-buffered EBSS (Earle's balanced salt solution) and some product codes are supplemented with 10 mg/l gentamicin. Sil-Select Stock can be used in combination with IUI, IVF and ICSI.

*For professional use only.*

### MATERIAL INCLUDED WITH THE KIT

Product code: SIS100

» 1x 100ml Sil-Select Stock

Product code: SIS500

» 1x 500ml Sil-Select Stock

Product code: SIS100\_32G

» 1x 100ml Sil-Select Stock with gentamicin

Product code: SIS500G

» 1x 500ml Sil-Select Stock with gentamicin

### MATERIAL NOT INCLUDED WITH THE KIT

- » 3cc syringes with 1 1/2" 21 g needle
- » Centrifuge (must be able to operate for up to 30 minutes at 400 g)
- » Incubator or water bath at 37°C (optional)
- » FertiCult Flushing Medium
- » LAF bench (ISO Class 5)

### PRODUCT SPECIFICATIONS

- » pH: 7.20-7.90 (Release criteria: 7.20-7.60)
- » Osmolality: 300-330 mOsm/kg
- » Density: 1,1150 – 1,1250 g/ml
- » Viscosity: < 1,75 cP
- » Endotoxin: < 0,5 EU/ml
- » Sterility: SAL 10<sup>-3</sup>
- » Sperm Survival test: ≥ 80% survival after 4 hours exposure of density selected spermatozoa to the test medium;
- » 75% survival after 24 hours exposure of density selected spermatozoa to the test medium
- » Not MEA tested
- » A certificate of analysis and MSDS are available upon request.

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### INSTRUCTIONS FOR PREPARATION OF GRADIENTS

We advise to prepare a dual gradient system (45% – 90% or 40% – 80%) starting from Sil-Select Stock. If preferred a multi-layer can be used as well (e.g. 45% – 70% – 90%). To prepare a 90% gradient mix 1 part FertiCult Flushing medium with 9 parts Sil-Select Stock. A 45% gradient is prepared by mixing 5.5 parts FertiCult Flushing medium with 4.5 parts Sil-Select Stock. Alternatively any HEPES-buffered EBSS-based medium can be used for the preparation of the gradients.

Note: Gradients should be prepared and repacked under sterile conditions (e.g. LAF bench ISO Class 5). For optimal results, prepare the gradients maximum 24 hours prior to use. Mix well after diluting the Sil-Select Stock solution. Follow the instructions as indicated for Sil-Select Plus in order to use the prepared density gradient media for sperm selection.

### STORAGE AND CONSERVATION

- » Store products without gentamicin between 2-25°C before first use, once opened store between 2-8°C.
- » Store products with gentamicin between 2-8°C.
- » Open and close bottles under aseptic conditions (e.g. LAF bench, ISO Class 5).
- » The products can be used safely up to 7 days after opening when sterile conditions are maintained and the products are stored at 2-8°C.
- » Do not freeze.
- » Do not use after expiry date
- » Content cannot be re-sterilized after opening
- » The products are stable after transport (max. 5 days) at elevated temperatures (≤ 37°C)

### WARNINGS AND PRECAUTIONS

All human, organic material should be considered potentially infectious. Handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.

Always work under strict hygienic conditions (e.g. LAF bench ISO Class 5) to avoid possible contamination, even when the products contain gentamicin. In case Sil-Select Stock medium contains gentamicin, appropriate precautions should be taken to ensure that the product is not sensitized to this antibiotic.

### TECHNICAL SUPPORT

FertiPro N.V.  
Industriepark Noord 32  
8730 Beernem / Belgium  
Tel +32 (0)50 79 18 05  
Fax +32 (0)50 79 17 99  
URL: www.fertipro.com  
E-mail: info@fertipro.com

0344

### STORAGE AND CONSERVATION

- » Store products without gentamicin between 2-25°C before first use, once opened store between 2-8°C.
- » Store products with gentamicin between 2-8°C.
- » Open and close bottles under aseptic conditions (e.g. LAF bench, ISO Class 5).
- » The products can be used safely up to 7 days after opening when sterile conditions are maintained and the products are stored at 2-8°C.
- » Do not freeze.
- » Do not use after expiry date
- » Content cannot be re-sterilized after opening
- » The products are stable after transport (max. 5 days) at elevated temperatures (≤ 37°C)

### WARNINGS AND PRECAUTIONS

All human, organic material should be considered potentially infectious. Handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.

Always work under strict hygienic conditions (e.g. LAF bench ISO Class 5) to avoid possible contamination, even when the products contain gentamicin. In case Sil-Select Stock medium contains gentamicin, appropriate precautions should be taken to ensure that the product is not sensitized to this antibiotic.

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