

To whom it concerns,

FertiPro N.V. updated the instructions for use of the Fructose Test. A summary of the changes is listed in the table below.

Instructions for use	Fructose Test
Previous version	FP09 I29 R01 B.10
New version	FP09 I29 R01 B.11
Changes	<ul style="list-style-type: none">- Small change in wording "absorbs" -> "can be measured"- Addition of reference to the website to download the MSDS and COA- Addition of 'freeze seminal plasma when it can not be tested at the same day'- Hyperlink to excel calculation sheet on website

Name: Laura Verspeelt
Function: Regulatory Affairs Associate
Date: 05/03/2019
Signature and stamp:

