

To whom it concerns,

FertiPro N.V. updated the instructions for use of the Citric Acid Test. A summary of the changes is listed in the table below.

| | |
|----------------------|---|
| Instructions for use | Citric Acid Test |
| Previous version | FP09 I37 R01 A.8 |
| New version | FP09 I37 R01 A.9 |
| Changes | <ul style="list-style-type: none">• Addition of blanc sample in the method• Addition of background correction in the calculation of results• Change in volume Reagent 1: 20ml -> 28ml• Change in volume Reagent 2: 10ml -> 11ml• Addition of chapter 'intended use' |

Name: Laura Verspeelt

Function: Regulatory Affairs Associate

Date: 29/04/2019

Signature and stamp:


**FertiPro nv**
Industriepark Noord 32
B-8730 Beernem / Belgium
VAT/TVA: BE 0448 767 035