

To whom it concerns,

FertiPro N.V. updated the instructions for use of VitalScreen. A summary of the changes is listed in the table below.

Instructions for use	VitalScreen
Previous version	FP09 I14 R01 C.6
New version	FP09 I14 R01 C.7
Changes	Addition of note: "In case nigrosin precipitates in the semen sample (rarely occurs in viscous samples), it is recommended to only work with Eosin (Reagent 1) and not to use Nigrosin (Reagent 2)"

Name: Laura Verspeelt  
Function: Regulatory Affairs Associate  
Date: 15/02/2019  
Signature and stamp:

