

To whom it concerns,

FertiPro N.V. updated the instructions for use of LeucoScreen Plus. A summary of the changes is listed in the table below.

Instructions for use	LeucoScreen Plus
Previous version	FP09 I111 R01 A.3
New version	FP09 I111 R01 A.4
Changes	<ul style="list-style-type: none">- Part of the introduction has been rewritten.- More information is provided on "40 days of analysis"- We do not advise to use a phase-contrast microscope- Reference is made to the website for the training program and demonstration video.- In case of very high concentration of round cells (i.e. above 20×10^6 per ml), it is strongly advised to dilute the sample in PBS or in FertiCult™ Flushing medium.

Name: Laura Verspeelt
Function: Regulatory Affairs Associate
Date: 3/02/2021
Signature and stamp:

 **Laura Verspeelt**
Regulatory Affairs
FertiPro N.V.

INVOICE ADDRESS

FertiPro nv
Industriepark Noord 32
B-8730 Beernem/Belgium

DELIVERY ADDRESS

FertiPro nv
Industriepark Noord 23
B-8730 Beernem/Belgium

tel. +32 (0)50 79 18 05
fax +32 (0)50 79 17 99
info@fertipro.com
www.fertipro.com

RPR Brugge
VAT BE 0448 767 035
IBAN BE 98 7374 4024 9293
BIC KREDBEBB