

To whom it concerns,

FertiPro NV updated the instructions for use of FertiCult IVF medium. A summary of the changes is listed in the table below.


Instructions for use	FertiCult IVF medium
Previous version	FP09 I07 R01 F.3
New version	FP09 I07 R01 F.5
Changes	<ul style="list-style-type: none"> <li>- Indications for use, intended user and target patient population are more clearly defined as separate items (regulatory requirement).</li> <li>- Optimization MEA protocol and addition of HSSA test as release criteria: change is reflected in instructions for use. <ul style="list-style-type: none"> <li>- Mouse Embryo Assay (% blastocysts after 96 hours, exposure to test medium: 24 hours): ≥80%</li> <li>- Human Sperm Survival Assay (% motility compared with control after 24 hours): ≥80%</li> </ul> </li> </ul>

Name: Eline De Reu

Function: Regulatory Affairs Associate

Date: 18/12/2025

Signature and stamp:



Eline De Reu  
Regulatory Affairs  
FertiPro N.V.

**COMPANY ADDRESS**

FertiPro nv  
Industriepark Noord 32  
B-8730 Beernem / Belgium

tel. +32 (0)50 79 18 05  
fax +32 (0)50 79 17 99  
info@fertipro.com  
www.fertipro.com

RPR Brugge  
VAT BE 0448 767 035  
IBAN BE 98 7374 4024 9293  
BIC KREDBEBB